




Part 1 Liter Vessel
 Part Number 02503
 Serial Number 96340 - 97276



		USP / Ph.Eur. Tolerances	
Dimensions	Vessel	Inside Ø a	100 mm 98 – 106 mm
		Height b	200 mm 160 – 210 mm
		With hemispherical bottom	
Material	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2024-01-11

Michael Schmidlin
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

Test Certificate