

# Test Certificate

Part: **Basket Holder with  
Stirring Shaft**

Number: 07558

Serial Number: N18002336 – N18002992

For SOTAX AT Dissolution Testing Systems



Dimensions:	Shaft Ø	9.4 to 10.1 mm	USP/Ph.Eur. Tolerances
	Cap Ø	25.0 mm	- -
	Cap size	5.1 mm	± 3.0 mm
	Clear opening	20.2 mm	± 0.5 mm
	Vent hole Ø	2.0 mm	± 1.0 mm
Material:	Stainless Steel		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management  
Bruno Balli  
2019-01-12

  
Signature