



Test Certificate

Part Protection Ring
Part Number 11148 and Set 11252 (7 pc)

For the Flange Protection of SOTAX Glass Vessels



		Tolerances		
Dimensions	Vessel	Outside Ø	125.6 mm	± 0.5 mm
		Inside Ø	124.0 mm	± 0.1 mm
		Inside Ø	40.0 mm	+ 0 / - 0.3 mm
		Height	5.95 mm	± 0.1 mm
		Max. temperature	60 °C	
Material	PP			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed.

Aesch, 2019-11-04

Bruno Balli
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland