




Test Certificate

Part 1 Liter Vessel Inactinic Brown
 Part Number 14143-02
 Serial Number P36 000016, P36 000019, P36 000020, P36 000023,
 P36 000025, P36 000027, P36 000030.



			USP / Ph.Eur. Tolerances	
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	98 – 106 mm
		Height b	168.2 mm	160 – 210 mm
		With hemispherical bottom		
Material	Borosilicate glass inactinic brown			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2021-04-30

4/30/2021 2:51:32 PM

SOTAX Quality Management
 EL 142/143 | Ground Floor | TTC Industrial area | MIDC Mahape | Navi Mumbai 400701 | India