




Part 1 Liter Vessel  
Part Number 14143-11  
Serial Number P35 001685 - P35 001845



		USP / Ph.Eur. Tolerances		
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	98 – 106 mm
		Height b	168.2 mm	160 – 210 mm
		With hemispherical bottom		
Material	Borosilicate glass			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2022-08-24

24-08-2022 15:20:22

Moreshwar Pulekar  
SOTAX Quality Management  
C-64 | First Floor | TTC Industrial area | Turbhe | Navi Mumbai 400705 | India

# Test Certificate