




Part 1 Liter Vessel
 Part Number 14322-01
 Serial Number 96644 - 96683



Test Certificate

		USP / Ph.Eur. Tolerances	
Dimensions	Vessel	Inside Ø a	98.5 to 101.5 mm
		Height b	200 mm
		With hemispherical bottom	160 – 210 mm
Material	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2023-10-19

Michael Schmidlin
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland