

Part Basket without Lock Mesh 10
 Part Number 15796-01
 Serial Number N3300088 – N33000113



		Tolerances		
Dimensions	Basket	Basket rim O.D. Ø	25.0 mm	± 3.0 mm
		Height	37.0 mm	± 3.0 mm
		Wire Ø	0.635 mm	--
		Wire openings	1.905 mm	--
		Screen opening height	27.0 mm	± 1.0 mm
		Screen O.D. Ø	22.2 mm	± 1.0 mm
		Clear opening Ø	20.2 mm	± 1.0 mm
Material	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2022-01-10

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland