




Part Plastic Sieve 10 mesh (3 Tubes)
 Part Number 23443-01



			Ph.Eur. Tolerances
Dimensions	A	Aperture size	1.8 to 2.2 mm
		Aperture shape	0.60 to 0.66 mm
Material	Plastic RTP 702 - HDPE		

Remark: Aperture shape square

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with Ph Eur. 2.9.1 Pharmacopeia requirements (at date of approval, see below).

Aesch, 2023-12-19

Michael Schmidlin
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

Test Certificate